

Mobile Watch System With GPS



Mobile Watch System with GPS

Getting To Know Your Mobile Watch System with GPS:

The Mobile Watch System with GPS is an affordable yet cutting edge MPERS solution that enables active customers to confidently engage with their world knowing that assistance is quickly available anywhere.

Runs on AT&T LTE Network

Nationwide mobile coverage supported by the widely adopted 4G LTE network.

All Day Battery Life

Best in class configuration for check-in and monitoring allows for up to 24 hours of battery life.

Proprietary Location Detection

The combination of GPS and WiFi provides the fastest and most accurate location service available.

Full Sensor Suite

Customers can monitor their heart rate and daily step count using our simple touch interface.

Vibration Feedback for Hearing & Visually Impaired

Don't let physical challenges hinder your customers from utilizing the full suite of services.

Coverage available in Canada and the USA.



GO ANYWHERE SAFELY

The Mobile Watch System with GPS is custom built for seniors who need a medical alert system to fit into their active lifestyles. It's an attractive and affordable option targeting a younger demographic than traditional PERS.

Specifications

The Mobile Watch System with GPS is ideal for seniors with an active lifestyle and individuals that spend a majority of their time on the go.



Cellular

- Supported Cellular Bands: LTE B2/B4/B5/B12
- Over The Air (OTA) configuration updates

Location Technology

- GPS (GPS & International GNSS)
- WiFi 802.11b/g/n
- LBS / Network based location backup

Physical & Environmental

- Dimensions: 47.5 mm x 16.9 mm
- Weight: 65 grams
- Screen: AMOLED True Color
- Water Resistance: IP67
- Operating Temperature: -20°C to +50°C
- Humidity: Up to 95% RH @ 50°C non-condensing

Electrical

- 5V DC 1A USB Charger
- 580mAh 3.7V Li-Po battery
- Battery Life up to 24 hours

Usability & Accessibility

- Intuitive user experience with voice prompts
- Accidental activations can be cancelled by user
- Haptic feedback with vibration motor
- Full sensor suite - heart rate, motion, vibration, and step counting
- Activity monitoring - users can see their daily step count
- User replaceable / customizable watchband

Charging Your Device

IMPORTANT: Before using your Mobile Watch System with GPS, it should be **fully charged**.



In Case of Emergency

If you need assistance, press and **HOLD** the Home Button until you hear the emergency alert message.

IMPORTANT: To cancel an accidental activation, press the **X** on the display prior to the countdown reaching zero.



This device is dependent on cellular coverage to work. Where cellular coverage is limited, you may encounter diminished device capabilities.

GETTING STARTED WITH YOUR MOBILE WATCH SYSTEM WITH GPS

Galaxy Medical Alert System Ltd.

Know Your Mobile Watch System with GPS

If the screen is blank, press the Home Button to display the Time & Date. When viewing the Time/Date, swipe right to access additional options.



Menu Options

Heart Rate

Press the watch face to measure your heart rate.



Pedometer

Press the watch face to track your daily steps and activity.



Turn Watch Off

Press the watch face to turn off your watch.



We recommend for safety purposes not to turn your watch off unless necessary.

Battery Warning

When the battery level is low, you will hear a low battery message and the battery level indicator will turn red. When this happens, you should charge your device as soon as possible.



Shower Safe

Your watch is water-resistant and can be worn while washing hands, watering plants and showering.

